

## CHOICE OF SERUM FOR EX VIVO EXPANSION DETERMINES IN VIVO CAR T CELL FUNCTION

Adapted from  
poster presentation at  
2019 3rd Annual Clinical Application of  
CAR T cells Conference  
at MSKCC.

..... Baylor College of Medicine, Houston, TX, USA  
..... INDIANAPOLIS, IN, USA

1102 Indiana Ave.  
Indianapolis, IN 46202

**SEXTON**  
BIOTECHNOLOGIES

**SEXTON**  
BIOTECHNOLOGIES

**T-LIVEN PR**

PATHOGEN-REDUCED HPL

EFFECTS ON CAR T CELLS

## T LIVEN IMPROVES TUMOR CELL KILLING FUNCTION OF CAR T CELLS


Adapted from poster presentation  
at 2018 ASGCT Annual Meeting,  
abstract # 468,  
Molecular Therapy, vol. 26,  
no. 5, pp. 220.

Center for Cell and Gene Therapy .....  
COOK REGENTEC .....




T-LIVEN PR SUPPORTS EXPANSION OF CAR T CELLS

- PSCA-CAR T (prostate stem cell antigen) cell generation was initiated in 10% fetal bovine serum (FBS) supplemented media.
- On Day 3 post transduction cells were split into 3 conditions: 10% FBS, 10% human AB serum (ABS) or 10% T-Liven PR and maintained with IL-2 (50U/mL) for 7 more days.
- Data are a mean ± S.E. of n=6 donors.




T-LIVEN PR IMPROVES MAINTENANCE OF LESS DIFFERENTIATED CAR T CELLS

- CAR T cell phenotype was analyzed measuring CCR7 expression by FACS on day 7 of expansion.
- The combined total percentage of T naïve (TN) and T central memory (TCM) cells for both CD4+ and CD8+ T cells are shown.
- Data are a mean ± S.E. of n=3 donors, \*:p<0.05 \*\*:p<0.01, \*\*\*:p<0.001.




T-LIVEN PR SUPPORTS IN VITRO CYTOTOXIC EFFECT OF CAR T CELLS

- The cytotoxicity of CAR T cells expanded for 7 days with different media supplements was tested against a PSCA positive cell line (Capan-1) using a <sup>51</sup>Cr-release assay.
- Data are a mean ± S.E. of n=4 donors, NT = nontransduced cells, CAR = PSCA transduced cells, E:T = Effector to Target cell ratio.




SOLID TUMOR: T-LIVEN PR IMPROVED IN VIVO CYTOTOXIC EFFECT OF CAR T CELLS

- Tumor cells injected SC (n=12/group)
- Animals treated with 1x10<sup>6</sup> CAR T cells grown in either AB Serum, FBS, or T-Liven
- Survival increased in T-Liven group




SOLID TUMOR: T-LIVEN PR SUPPORTS IN VIVO CYTOTOXIC EFFECT OF CAR T CELLS IN TUMOR RE-CHALLENGE MODEL

- Experimental design, re-challenge model.
- T cell bioluminescence on day of re-challenge.
- Re-challenge tumor measurement. Arrow indicates tumor re-challenge



LIQUID TUMOR: T-LIVEN PR IMPROVED IN VIVO CYTOTOXIC EFFECT OF CAR T CELLS

- Treated mice re-challenged with new tumor injection
- T-Liven group had higher residual CAR T cells at re-challenge
- T-Liven group maintained better tumor killing effect after re-challenge



CONCLUSION

THIS STUDY SHOWS THAT EXPOSURE OF CAR T CELLS TO T-LIVEN PR DURING EX VIVO CULTURE LEADS TO ENHANCED PERSISTENCE, BETTER IN VIVO ANTI-TUMOR ACTIVITY, AND SURVIVAL.